

May Livzon Bring The Reviving

Hope to The World

***Frontier Standard Traceability
(WHO International Standard)
Excellent Performance***

Evaluation solutions for
COVID-19 vaccines



PRODUCT OVERVIEW

The product can be used to evaluate the immune effect after vaccination or whether neutralizing antibodies are produced in human body after infection with 2019-nCoV.

PRODUCT APPLICATION

- Vaccine effectiveness evaluation
- Protective immunity indication
- Revaccination assessment
- Long-term protective immunity after recovery

KEY FEATURES

• Clinical performance

Compared with reference test by testing 125 positive and 147 negative specimens for SARS-CoV-2 neutralizing antibody.

Method	Reference Test		Total
	Positive	Negative	
2019-nCoV Neutralizing Antibodies Test (Fluorescence Immunochromatographic Assay)	123	1	124
	2	146	148
Total	125	147	272


Positive percent agreement: $123/(123+2) \times 100\% = 98.4\%$ [95% CI: 94.4%–99.9%]
Negative percent agreement: $146/(146+1) \times 100\% = 99.3\%$ [95% CI: 95.2%–99.9%]
Overall percent agreement: $123/(146+146) \times 100\% = 98.9\%$ [95% CI: 95.8%–99.6%]

• Limit of detection*

Not more than 5 IU/mL

*Traceability: First WHO International Standard for anti-SARS-CoV-2 Immunoglobulin (human) NIBSC code: 20/136

• Linearity



ORDER INFORMATION

Product Name	Cat. No.	Shelflife	Specimen	Storage Conditions	Packaging
2019-nCoV Neutralizing Antibodies Test	G01802020	12 months	Serum/Plasma/ Venous whole blood	2~30°C	20 tests/kit