
COVID-19 ANTIGEN RAPID TEST (Colloidal Gold)

P/N: FGCOVG300 (Sample type: Nasopharyngeal Swab)

P/N: FGCOVG1300 (Sample type: Anterior Nasal Swab)

P/N: FGCOVG1600 (Sample type: Oropharyngeal Swab)

JOINSTAR BIOMEDICAL TECHNOLOGY CO., LTD.

CONTACT:

Mobile: +86-13957167265 (Stanley) E-mail: info@joinstar.cn Tel: +86-571-89023160 Fax: +86-571-89028228

Web: <http://en.joinstar.cn>

 @JoinstarBioTech

 @JoinstarBioTech

 www.linkedin.com/company/joinstar

COVID-19 ANTIGEN RAPID TEST (Colloidal Gold)

P/N: FGCOVG300 (Sample type: Nasopharyngeal Swab)

- **Sample type:** Nasopharyngeal Swab
- **High sensitivity & specificity.** Can be used for early screening. Facilitates patient treatment decisions quickly;
- **Suitable for multiple scenarios:** medical institution testing; screening before resumption of work and school, continuous monitoring, etc.
- **One-step method,** easy to operate, reducing missed or false inspections caused by operator errors;
- **All necessary reagent provided & No equipment required;**
- **Time saving procedures,** results are available in 10 minutes;
- **Storage temperature:** 2~30 °C . No cold-chain transportation needed;
- **Specification:** 25 tests/box, 1 test/box;
- **Diverse cooperation modes:** OEM/ODM accepted.

COVID-19 ANTIGEN RAPID TEST(Colloidal Gold)

P/N: FGCOVG1300 (Sample type: Anterior Nasal Swab)


P/N: FGCOVG1600 (Sample type: Oropharyngeal Swab)

- **Sample type:** Anterior Nasal Swab&Oropharyngeal swab
- **Non-invasive, painless sample collection**
- **High sensitivity & specificity.** Can be used for early screening. Facilitates patient treatment decisions quickly;
- **Suitable for multiple scenarios:** medical institution testing; screening before resumption of work and school, continuous monitoring, etc.
- **One-step method,** easy to operate, reducing missed or false inspections caused by operator errors;
- **All necessary reagent provided & No equipment required;**
- **Time saving procedures,** results are available in 15 minutes;
- **Storage temperature:** 2~30 °C . No cold-chain transportation needed;
- **Specification:** 25 tests/box, 1 test/box;
- Diverse cooperation modes: OEM/ODM accepted.


Performance Characteristics


Sample type: Nasopharyngeal Swab

95% CI: 91.64%~99.39%
Sensitivity


95% CI: 97.96%-99.78%
Specificity


95%CI: 97.62%-99.53%
Accuracy


Sample type: Anterior Nasal Swab&Oropharyngeal swab

95% CI:(90.26%-98.92%)
Sensitivity



95%CI: (97.96%-99.78%)
Specificity

95%CI: 97.40%-99.42%
Accuracy

Nasopharyngeal Swab Sample collection

P/N: FGCOVG300 (Sample type: NasopharyngealSwab)



1. Carefully insert the swab into the nostril of the patient, reaching the surface of posterior nasopharynx. that presents the most secretion under visual inspection

2. Swab over the surface of the posterior nasopharynx. Rotate the swab 3~4 times

3. Withdraw the swab from the nasal cavity

Anterior Nasal Swab Sample collection

P/N: FGCOVG1300 (Sample type: Anterior Nasal Swab)


1. Tilt the patient' s head back 70 degrees.

2. While gently rotating the swab insert swab about 2~2.5 cm into nostril



3. Rotate the swab five times against the nasal wall then slowly remove from the nostril. Using the same swab repeat the collection procedure with the second nostril.

Oropharyngeal Swab Sample collection

P/N: FGCOVG1600 (Sample type: Sample type: Oropharyngeal Swab)

1. Tilt the patient's head slightly and make "ahhh" sounds, exposing the pharyngeal tonsils on both sides

2. Hold the swab and wipe the pharyngeal tonsils on both sides of the patient with moderate force back and forth for **at least 3 times**

Test Procedure


for COVID-19 Antigen Rapid Test (Colloidal Gold)

P/N: FGCOVG300 (Sample type: Nasopharyngeal Swab)

P/N: FGCOVG1300 (Sample type: Anterior Nasal Swab)

P/N: FGCOVG1600 (Sample type: Oropharyngeal Swab)

Please read the instructions carefully and allow the test device and specimens to equilibrate to temperature (15°C-30°C) prior to testing.

1. Open the cap of the extraction tube


2. put the anterior nasal swab with the sample into the extraction tube and rotate it 10 times

3. Squeeze the wall of the extraction tube by hand


Test Procedure

for COVID-19 Antigen Rapid Test (Colloidal Gold)


Please read the instructions carefully and allow the test device and specimens to equilibrate to temperature (15°C-30°C) prior to testing.

4. Take out the swab



5. Cover the cap of the extraction tube, set aside.

OR

6. Cut off the protrusion of the collection tube **OR** unscrew the little transparent cap

Positive

Negative

Invalid

7. Take the test cassette from the packaging bag, place it on a table

8. Add 2 drops of the sample into the sample hole vertically

9. Read the result after 15 minutes. If left unread for 20 minutes or more the results are invalid and a repeat test is recommended.


CERTIFICATE

Business license

Medical device production license

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CE

DECLARATION OF CONFORMITY

Manufacturer: Joinstar Biomedical Technology Co.,Ltd.

Address: 10th Floor, Administration Building, NO.519, XingGuo RD., Yuhang Economic and Technological Development Zone, Hangzhou, Zhejiang, China, 311188

EC Representative's Name: Lotus NL B.V.

EC Representative's Address: Koningin Julianaplein 10, 1e Verd, 2595AA, The Hague, Netherlands.

Declares, that the product

Product Name and Model:

COVID-19 Antigen Rapid Test (Colloidal Gold)

1Test/Kit, 20Tests/Kit, 25Tests/Kit

as described above are in conformity with the requirements as defined in IVDD98/79/EC Annex III.

Additional information:

Conformity assessment route: Directive 98/79/EC, Annex III

Classification: List Others

I, the undersigned, hereby declare that the medical devices specified above conform with the Directive 98/79/EC on in vitro diagnostic medical devices and pertinent essential requirements.

Date Signed:

2020.11.02

Zhong Wang
Management Representative
Joinstar Biomedical Technology Co., Ltd.

Joinstar Biomedical Technology Co., Ltd.

Version: 0.0

CE Declaration of Conformity

COVID-19 In Vitro Diagnostic Devices and Test Methods Database

WHO-FIND List

Registration in Italy

- [Jingpo Microing Biomedical Technology Co., Ltd](#) Novel coronavirus 2019-nCoV nucleic acid detection kit (RT-PCR method) [\(PRK\)](#)
- [Jingpo Micro Biom Biotechnology Co., Ltd](#) Fluorescent RAA Detection for 2019-nCoV [\(PRK\)](#) [Contact](#)
- [Jingpo Biomedical Biotechnology \(Nanjing\) Co., Ltd](#) SARS-CoV-2 (COVID-19) IgM/IgG Antibody Fast Detection Kit (Colloidal Gold) [\(PRK\)](#) [Contact](#)
- [Jingpo Yifu Medical Apparatus](#) Detection Kit for Novel Coronavirus 2019-nCoV RNA (RT-PCR) Fluorescence Probing [\(PRK\)](#) [Contact](#)
- [Jinhuan Medical Instrument Co., Ltd](#) (COVID-19) IgM/IgG Antibody Fast Detection Kit (Colloidal Gold) [\(CE-IVD\)](#) [Contact](#)
- [JIN Winlays](#) PreTest Covid-19 RT-qPCR kit [\(Singapore HSA - CE-IVD\)](#) [Contact](#)
- [Joinstar Biomedical Technology Co., Ltd](#) COVID-19 Antigen Test (Latex) (CE-IVD) [Contact](#)
- [Joinstar Biomedical Technology Co., Ltd](#) COVID-19 Antigen Rapid Test (Colloidal Gold) (CE-IVD) [Contact](#)
- [Joinstar Biomedical Technology Co., Ltd](#) SARS-CoV-2 IgM/IgG Antibody Test (Colloidal Gold) (CE-IVD) [Contact](#)
- [Joinstar Biomedical Technology Co., Ltd](#) SARS-CoV-2/MERS-CoV/ Influenza A&B Antigen Rapid Test (CE-IVD) [Contact](#)
- [JQIEMIG \(Tianjin\) Biotechnology Co., Ltd](#) COVID-19 (SARS-CoV-2) Antigen Rapid Test Kit (Colloidal Gold) [\(CE-IVD\)](#) [Contact](#)
- [JQIEMIG \(Tianjin\) Biotechnology Co., Ltd](#) COVID-19 (SARS-CoV-2) IgM/IgG Rapid Test Kit (Colloidal Gold) [\(CE-IVD\)](#) [Contact](#)
- [JQIEMIG \(Tianjin\) Biotechnology Co., Ltd](#) COVID-19 IgM/IgG Antibody Test Kit (Latex) [\(CE-IVD\)](#) [Contact](#)
- [Kephene Diagnostics](#) KID-Rapid SARS-CoV-2 Antigen Test [\(in development\)](#) [Contact](#)
- [Kephene Diagnostics](#) KID-COVID-19 IgM/IgG Rapid Detection Test Kit [\(in development\)](#) [Contact](#)
- [Kephene Diagnostics](#) KID-SARS-CoV-2 IgM/IgG ELISA [\(in development\)](#) [Contact](#)
- [KIT Medical Co., Ltd](#) RND-COVID-19 Detection Kit and RND-COVID-19 Triple Detection Kit [\(CE-IVD\)](#) [Contact](#)
- [Kitch Biotechnology \(Nanjing\) Co., Ltd](#) SARS-CoV-2 Antigen/ Latex/Flow Assay [\(SARHA-UK\)](#) [Contact](#)

The database contains publicly available In Vitro Diagnostic Medical Devices for COVID-19 and it is updated periodically. Please note that additional performance (as retrieved from manufacturers web pages) is provided only for devices commercially available with CE-IVD mark.

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Have passed the performance verification of Paul Ehrlich Institute in Germany

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